



Single-Use Bag Totes

By G & G Technologies

Flexible design for custom and standard applications of buffer and media preparation, product formulation, and mixing of process fluids.

- 304, 316L stainless steel or polypropylene
- Jacketed and non-jacketed options for heating and cooling options
- Load cell options
- Cylindrical, quadrangular and multi-faceted

Single-Use Bags are becoming the alternative to expensive methods such as CIP, raising process efficiency and reducing batch-to-batch cross contamination with rapid turn-over time.



Standard Alternative Features include:

Weight Monitoring
Temperature Monitoring and Control
pH/ DO/ Conductivity Monitoring
Agitation Control

Other optional accessories are:

Lexan Top Covers
View/ Access Port Windows
Probe Managers
Test Tube holders
Custom Instrument Mounting

Standard Sizes: 20L, 50L, 100L, 200L, 500L, and 1000L

For more information, Contact G&G Technologies at
info@ggtechnologies.com or call us at 401-295-4000