



G&G TECHNOLOGIES, INC.

6 Grandview St.
Coventry, RI 02816
401-295-4000

www.ggtechnologies.com


Tubing Management Racks and Stands

Tubing Management

Whether you are building a new process or outfitting an existing one, tubing management is an important consideration when it comes to safety, cleanability, and product recovery. Keep your single-use flow-paths organized, elevated, and drainable, with our wide array of tubing management solutions.

Specific Features

- Pivot Capability
- Modular Design
- Telescoping Assembly
- Pitched runs for draining
- Mobile and Stationary Arrangements
- Stainless Steel Construction
- Variety of Styles



Advantages of G&G Technologies Tubing Management Systems

At G&G Technologies, we recognize the critical importance of durability in laboratory and production environments. Our tubing management components are crafted from high-quality materials, including 304 and 316 stainless steel, Delrin, and HDPE. These materials ensure longevity, corrosion resistance, and the ability to withstand the rigorous demands of production in life-sciences settings.

When it comes to tubing management, customization is key. G&G Technologies is dedicated to providing tailored solutions that address your specific tubing management challenges. We collaborate closely with you to understand your unique requirements and deliver systems that are flexible, efficient, and add value to your process.



Eliminate tubing management bottlenecks in your facility. Contact G&G Technologies today, and let our experts create a customized Tubing Management System that optimizes your processes and elevates your pharmaceutical production. Efficiency, flexibility, and precision are within reach with G&G Technologies.



Mobile Tubing Management Stands

The G&G Mobile Tubing Management Stands are a flexible solution to improve tubing management in your facility. These stands are designed for mobility, allowing you to easily relocate them within your facility. Their compact and mobile design makes them convenient to reconfigure your production line and address workflow challenges.



















These stands serve as a versatile platform for tubing management guides, creating a dynamic routing system tailored to your operational needs. The stainless-steel rails provide a stable foundation, and central hinges allow pivotal movement, enhancing adaptability in navigating tight spaces or accommodating changes in tubing direction. The telescoping assembly takes mobility further, allowing you to route around obstacles. Whether routing tubing over or under equipment, the telescoping feature ensures a downward sloping profile, optimizing fluid transfer.





Mobile Tubing Management Stands Product Array
















Explore our robust product array of tubing management stands, rail systems, and guides. Our offerings include a variety of form factors, in a variety of sizes, providing a stable and mobile base for the rail systems. H-Carts are available in a variety of sizes to accommodate tubing runs that may be sparsely or densely populated. Delta Carts are designed to nest together, allowing for compact storage when not in use. Complementing these are our hinged assemblies and junctions, which fasten the rail systems in place. Our product offering features a diverse array of tube guides and anchors that mount onto the rail systems, allowing secure and precise placement throughout your facility. With this diverse range of products, G&G Technologies provides a complete tubing management system, offering efficiency and customization tailored to the specific demands of your single-use facility.

Carts		Cart Accessories		Rail-Mount Tube Guides - Examples		
					5-Wide	10-Wide
Jumbo H Carts		Rail Junction		Slip-On		
	F491479		F491371			
H Carts		Rail Junction Tee		Slip-On Locking		
	F491481		F491486			
Mini H Carts		Rail Junction X		Clamp-On		
	F491561		F492084			
Delta Carts		Top guide		Clamp-On Locking		
	F491586					
		Clamp-On Guide		Custom configurations are available for each style		
		T-Loc Panels				



Stationary Tubing Management

Introducing our Stationary Tubing Management System, designed for permanent installation in your facility. These stands are fastened securely in place, ensuring a sturdy setup for tubing management. The bollards secure rail systems for tubing management as well as conduits for electrical routing. Additional accessories include holders for transmitters and valves, for flexible placement in your facility.

Bollards		Rail-Mount Tube Guides - Examples			Accessories	
			5-Wide	10-Wide		
Bollard Single Run		Slip-On			Transmitter Holders	
Bollard Double Run		Slip-On Locking			Valve Holders	
Bollard Corner		Clamping			Sensor Holders	
Bollard Tee Junction		Clamping Locking			Filter Holders	